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SUBJECT

Cleaning and sterilization of non-sterile medical devices of Spierings Orthopaedics. If specific medical devices are excluded it will be mentioned in the instructions for use of the medical device.

WARNING

- We refuse any liability regarding the quality of the instruments if they have been used on patients with Creutzfeldt-Jakob disease or HIV infection disease.
- Soiled medical devices will not be cleaned at the position of the clamping point with the fixation strips when the instruments are placed in the sterilization tray during cleaning.
- If the conditions do not allow medical devices to be cleaned outside the fixation, the user must ensure that the medical devices are clean at the fixation points before they are placed in the sterilization process.
- Never use metal brushes or steel wool for cleaning.
- Do not exceed a temperature of 140 °C during reprocessing.
- Prolonged immersion of the device in a sodium chloride solution (NaCl) may cause pitting or stress corrosion.
- Excessive pressure, twisting or inappropriate use of the medical device can damage it.
- Never use damaged medical devices. Repairs must only be made by the manufacturer. Repairs made by anyone else, void the guarantee and CE marking.
- The use of paraffin oils on plastic (parts of) devices may cause absorption of paraffin into the device resulting in swollen parts.
- Refurbishing of the device is only allowed by the manufacturer.

GENERAL

To minimize the risk for the patient, these instructions should be carefully followed. The correct application, the use of suitable equipment, the right procedures and the execution by qualified and properly trained personnel, is the full responsibility of the user. These instructions for use have been validated in accordance with ISO 17664. It remains the responsibility of the facility to ensure that the procedure is performed using equipment, materials and personnel at a designated area, and achieves the desired requirements. This includes validation and routine monitoring of the process. Likewise, any deviation by the processor from the described procedures should be evaluated for effectiveness and any potential adverse consequences.

- Reprocessing must start immediately after use.

Procedure

1. Pre-cleaning
2. Automated cleaning and disinfection
3. Inspection and testing
4. Packaging
5. Sterilization
6. Storage

1. PRE-CLEANING

- Pre-cleaning must start immediately after use to prevent fixation of proteins and other contaminants.
- If applicable medical devices should be disassembled according to the Surgical Technique of Spierings Orthopaedics BV.
- Detergents without protein fixing effect should

be used. Follow the detergents manufacturer's instructions.

- Only use ultrasonic devices which are intended for use with medical devices.
1. Remove surface debris under running water.
 2. Apply the detergent to the complete surface of the medical device, using a clean brush. Medical devices with joints have to be cleaned in an open and closed position. Use appropriate brushes for cleaning the inside of narrow lumen and blind holes.
 3. Hold the medical device under running water and rinse thoroughly, use a syringe for narrow lumen. This should be repeated several times.
- In case of narrow lumen or blind holes subsequent ultrasonic cleaning is required.

Ultrasonic cleaning

4. Preheat the ultrasonic bath to the temperature recommended by the manufacturer of the ultrasonic bath or detergent (if used).
5. Completely immerse the medical device in the cleaning liquid. Narrow lumen and blind holes must be filled with the liquid. Water should be used as a cleaning liquid, a detergent may also be used.
6. Clean the medical device in the ultrasonic bath.
7. Flush the outer surface as well as the narrow lumen and blind holes using a syringe with water. This should be repeated several times.

2. AUTOMATED CLEANING AND DISINFECTION

- Cleaning agent: Neodisher MediClean (forte) of Chemische Fabrik Dr. Weigert & CO.KG, Hamburg, Germany or similar, 2-10 ml/l (0.2-1.0%).
- A washer-disinfector meeting the requirements of ISO 15883 or equivalent standard is required.

1. Medical devices with joints must be placed with the hinges open. Narrow lumen medical devices and other medical devices with hard to get or blind holes must be positioned in a way that water is able to run thoroughly through the medical device and assist drainage.
2. The cycle must include: pre-rinse cycle with cold water, wash at 40-60 °C for at least 10 minutes, rinse with cold purified water (e.g. RO water).
3. Disinfect at 90 °C for at least 5 minutes.
4. Dry using hot air of 99 °C for at least 35 minutes.
5. Check the medical device, especially narrow lumen and blind holes, to insure proper cleaning. If not, repeat the cleaning and disinfection process or clean manual/ultrasonic before repeating the automated cleaning and disinfection.

3. MAINTENANCE AND INSPECTION

Verify that the instrument's form and function have not deteriorated due to repeated use and reprocessing. The medical devices should be checked to assume that all parts are moving correctly. It is allowed to use a droplet of high quality oil on hinges. This oil must be physiologically safe (e.g. paraffin oil according to DAB 8 of Ph. Eur. or USP XX). The use of paraffin oils on plastic (parts of) devices may cause absorption of paraffin into the device resulting in swollen parts. A visual check should be done on the medical device to see if there is any damage. Each medical device that is part of a large set, should be checked with all the other parts.

4. PACKAGING

Individual:

A standardized packing material (e.g. EN 868) can be used. Avoid any tension on the packaging material when it is sealed.

Sets:

All the medical devices must be placed in the correct location of a customized sterilization tray.

5. STERILIZATION

Sterilize using moist heat in accordance with EN ISO 17665 at 134 °C with an exposure time of at least 3 minutes (the medical devices can withstand sterilization up to 18 minutes). If other sterilization methods are used, Spierings Orthopaedics BV should be contacted and this method must be validated by the sterilization institution. The device should be inspected during sterilization for indications of deterioration which could compromise the fit and function of the device. These indications include: broken, deformed, does not fit in counterpart, corrosion, loss of sharpness, burrs, and wear. In case any indications of deterioration are present, the device must be taken out of service. After sterilization check all packaging for any damage and dryness and the sterility indicators for the right colour change.

6. STORAGE

Store sterilized medical devices in closed cabinets or sealed cabinets, protected against dust, humidity and temperature change. A maximum shelf life before use should be defined by each health care facility.

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Spierings Orthopaedics B.V.